


OncoMonitor®MRD is a blood test designed for Longitudinal Cancer Monitoring and Minimal Residual Disease Detection, by measuring CTC & ctDNA.



OncoMonitor® TRM is a liquid biopsy-based solution designed for monitoring treatment response, emergence of resistance and tracking disease progression.

## OncoMonitor®MRD Benefits



**Minimal Residual Disease Detection**

OncoMonitor's ability to identify trace amounts of ctDNA and CTCs allows for early detection of residual cancer cells post-treatment, enabling timely intervention and improved patient outcomes.



**Prognostic Value**

CTC and ctDNA levels have shown prognostic significance in various cancers. The test results can help predict disease progression and overall survival, aiding in informed decision-making for both patients and healthcare providers.

### ctDNA epigenomics validation results

Sensitivity	90.5%
PPV	95%
Specificity	97.2%
NPV	94.6%
Concordance	94.7%
Limit of Detection (LOD)	0.05

Analysis done on 57 samples which includes Seraseq® reference standards and healthy individual samples

## OncoMonitor® TRM Benefits



**Longitudinal Monitoring**

By regularly assessing CTC numbers and ctDNA load, oncologists can gain insights into the dynamics of the disease, its response to treatment, and potential disease relapse.




**Treatment Efficacy Assessment**

Changes in CTC counts and ctDNA levels can indicate whether the current treatment approach is successful or if modifications are necessary to enhance its impact.

## For which cancer phase is OncoMonitor® most applicable?

Applicability Across Cancer Journey: Tracking Disease Progression and Relapse Detection Through Longitudinal Blood-Based Monitoring



# Comparison of OncoMonitor® MRD and OncoMonitor® TRM


	OncoMonitor® MRD	OncoMonitor® TRM
Eligibility	Early-stage post-curative intent surgery	Advanced metastatic cancers
Features	Cellulomics and Epigenomics	Cellulomics, Epigenomics, and Genomics (mutation profiling)
Clinical Decisions	1. Escalation or de-escalation in adjuvant setting 2. Surveillance of tumor burden 3. MRD assessment	1. Early course correction 2. Predicting therapy response 3. Disease progression prediction
Test Type	Prognostic/Predictive test	Predictive test
Indications	Colon, Rectal, Breast, Bladder, Prostate, and Ovary cancers	All Solid Tumors



## Test Performance Characteristics:


- **Genomic Alterations Coverage:** Detects SNVs, INDELs, CNVs, and selected fusions (by DNA) in targeted genes (exclusive to TRM).
- **AI-Powered Reporting:** Enabled via iCARE™.
- **Sequencing Depth:** 10,000x mean depth.
- **Gene Coverage:** 99% mean coverage across all included genes.
- **Fusion Gene Analysis:** Covers 30 fusion genes with 40 corresponding partner genes, including DNA-based fusions.
- **Variant Calling Databases:** Utilizes approximately 20 databases.
- **MSI Status:** Assessed by NGS.
- **PD-L1 Expression Analysis:** Evaluated via immunofluorescence on OncoDiscover-CTC.
- **Indications:** Solid tumors.
- **CTC Detection:** Enabled through DCGI-approved OncoDiscover® technology.
- **DNA Methylation Profiling:** Covers over 3,440 CpG sites.
- **Sample Type & Volume:** Peripheral blood | 3 tubes (10ml PAXgene in 2 tubes, 10ml K2 EDTA in 1 tube).

## Key Actionable Genes Include in OncoMonitor TRM:




**Lung**

AKT1, ALK, BRAF, ERBB2 (Her2), EGFR, FGFR 1/2/3, HRAS, KRAS, MAP2K1, MET, NRAS, NTRK 1/2/3, PIK3CA, PTEN, RET, ROS1, TP53




**Breast**

AKT1, AR, BRCA1, BRCA2, ERBB2 (Her2), ESR1, FGFR 1/2/3, NTRK 1/2/3, PALB2, PIK3CA, PTEN, TP53




**Melanoma**

BRAF, CTNNB1, GNA1, GNAQ, HRAS, KIT, KRAS, MAP2K1, NF1, NRAS, NTRK 1/2/3, PDGFRA, PIK3CA, PTEN, TP53




**Bladder**

MLH1, MSH2, MSH6, NTRK 1/2/3, PMS2, TP53, TSC1




**Colorectal**

AKT1, BRAF, HRAS, KRAS, MET, MLH1, MSH2, MSH6, NRAS, NTRK 1/2/3, PIK3CA, PMS2, PTEN, SMAD4, TP53



**Ovary**

BRAF, BRCA1, BRCA2, HRAS, KRAS, NRAS, NTRK 1/2/3, PALB2, PDGFRA, TP53



**Gastrointestinal**

BRAF, KIT, KRAS, MET, MLH1, MSH2, MSH6, NTRK 1/2/3, PDGFRA, PMS2, SMAD4, TP53

Pioneering Precision Monitoring  
Navigating Cancer’s Ebb and Flow

Longitudinal Monitoring

Assess disease dynamics over time

Minimal Residual Disease detection

Uncover residual cancer cells post-treatment

AI powered Real-time results

Swift and accurate AI powered comprehensive reporting

